

20ml to 60ml Vials

- Volume from 20ml to 60ml
- PTFE septa can resist organic solvent
- Caps are made of high quality polypropylene



Specifications

Model No.	Size	Description
V20A22ET	27x58mm	20ml screw-thread vial, clear glass
V20B22ET		20ml screw-thread vial, amber glass
SC22WPEET	Φ22mm	White closed top PP cap + 22mm white PE septa
V20A	27x58mm	20ml screw-thread vial, clear glass
V20B		20ml screw-thread vial, amber glass
V30A	27x75mm	30ml screw-thread vial, clear glass
V40A	27x95mm	40ml screw-thread vial, clear glass
V40B		40ml screw-thread vial, amber glass
V60A	27x140mm	60ml screw-thread vial, clear glass
V60B		60ml screw-thread vial, amber glass
SC22A22A	Φ22x3mm	nature PTFE/nature Silicone septa 22*3mm, 22mm white screw PP cap, 15mm centre hole
SC22A22B	Φ22x1.5mm	nature PTFE/nature Silicone septa 22*1.5mm, 22mm white screw PP cap, closed-top
SC22A22C	Φ22x3mm	nature PTFE/nature Silicone septa 22*3mm, 22mm black screw PP cap, 15mm centre hole
SC22A22D	Φ22x1.5mm	nature PTFE/nature Silicone septa 22*1.5mm, 22mm black screw polypropylene cap, closed-top
S22A-I	Φ22x3mm	nature PTFE/nature Silicone septa 22*3mm
S22A-II	Φ22x1.5mm	nature PTFE/nature Silicone septa 22*1.5mm