

# 2ml Vials & Caps (8-425)

- High quality
- Unique thread design
- USP TYPE I Grade
- OEM Accepted



### Specifications

Model No.	Size	Description
V8A	11.6x32mm	2ml small opening short screw-thread vial,clear
V8B		2ml small opening short screw-thread vial with write-on spot,clear
V8C		2ml small opening short screw-thread vial,amber
V8D		2ml small opening short screw-thread vial with write-on spot,amber
SC8A8A	Φ8mm	white PTFE/red Silicone septa, black screw PP cap,5.5mm centre hole
SC8AA8A		pre-slit white PTFE/red Silicone septa, 8mm black screw PP cap,5.5mm centre hole
SC8F8F		Nature PTFE/white Silicone spata, 8mm blace screw PP cap, 5.5mm centre hole
S8A		white PTFE/red Silicone septa
S8F		Nature PTFE/white Silicone speta

100pcs/pk