

# 9mm Micro Vials

## Glass



### Specifications

Model No.	Size	Description
VL9A	9x11.6mm	300ul thread fusion type vials, clear, borosilicate.
VL9B	9x11.6mm	300ul thread fusion type vials, clear, write-on spot, borosilicate.
VL9C	9x11.6mm	300ul thread fusion type vials, amber, borosilicate.
VL9D	9x11.6mm	300ul thread fusion type vials, amber, write-on spot, borosilicate.

100pcs/pk

## Plastic



### Specifications

Model No.	Size	Description
VP9A	9x11.6mm	0.3ml Short thread micro injection sample vials, clear, Polypropylene
VP9C	9x11.6mm	0.3ml Short thread micro injection sample vials, amber, Polypropylene

100pcs/pk