

# Anti-bacterial Filter

## 50mm

### Product Description

#### PASS HEPA TEST

Accepted ID 7-13mm tubing

High flow rates at low differential pressures

Rugged PC housing withstands high temperatures, pressures and aggressive solvents

High particulate-loading capability to maintain flow rates over extended periods

Manufactured in clean room under controlled Conditions

Autoclave



### Application:

1. In-line gas filter
2. Deliver Sterile air
3. Smoke excluded
4. Benchtop Fermentors
5. Vacuum Pump
6. liquid filtration, media membrane including PVDF、Hydrophilic PTFE、Nylon、MCE、PES for customer choose

