

### A Test Report Demonstrating Nuclease Content of Plastic Products

- Requester: Finetech Research and Innovation Corporation
- Analyzer: Medical Analysis Lab,
   Department of Medical Laboratory Science and Biotechnology,
   Yuanpei University of Medical Technology
- Goals of the tests:

By using PCR methods, our goals of the tests were to prove that centrifuge tubes manufactured by Finetech are both DNase-free and RNase-free.



#### Outline

- Goals of the tests
- A test method analyzing DNase content of centrifuge tubes
- Test results 50 mL centrifuge tubes
- Test results 15 mL centrifuge tubes
- A test method analyzing RNase content of centrifuge tubes
- Test results 50 mL centrifuge tubes
- Test results 15 mL centrifuge tubes
- Conclusions

## Test results - 50 mL centrifuge tubes




Lanes 50 mL tubes Supplier A(P) 1 Supplier A(N) 2 **Finetech** 3 Finetech 4 5 Supplier B 6 Supplier C 7 Supplier D

Test results:

50 mL centrifuge tubes of Finetech are DNase-free

# Test results -15 mL centrifuge tubes



15 mL tubes

Supplier A(P)

Supplier A(N)

**Finetech** 

Finetech

Supplier B

Supplier C

Supplier D

Lanes

1

2

3

4

5


6

7

Test results:

15 mL centrifuge tubes of Finetech are DNase-free

## Test results - 50 mL centrifuge tubes




| Lanes | 50 mL tubes   |  |
|-------|---------------|--|
| 1     | Supplier A(N) |  |
| 2     | Supplier A(P) |  |
| 3     | Finetech      |  |
| 4     | Finetech      |  |
| 5     | Supplier B    |  |
| 6     | Supplier C    |  |
| 7     | Supplier D    |  |

Test results:

50 mL centrifuge tubes of Finetech are RNase-free

### Test results - 15 mL centrifuge tubes



| Lanes | 15 mL tubes   |  |
|-------|---------------|--|
| 1     | Supplier A(N) |  |
| 2     | Supplier A(P) |  |
| 3     | Finetech      |  |
| 4     | Finetech      |  |
| 5     | Supplier B    |  |
| 6     | Supplier C    |  |

Test results:

15 mL centrifuge tubes of Finetech are RNase-free



#### Conclusions

#### Results of the tests:

• According to these test results, it was proven that centrifuge tubes manufactured by Finetech are both DNase-free and RNase-free.

• The quality of tubes manufactured by Finetech was assured to be as good as that of tubes produced by some well-known worldwide brands which were used as controls.

| Brand    | DNase | RNase |
|----------|-------|-------|
| Finetech | _     | _     |