

Micro Insert & Shell Vials



Specifications

Model No.	Size	Description
IPA150	29x5mm	150ul Insert with mandrel interior and polymer feet, suits for ND8 Vials
IPB250	29x5.7mm	250ul Insert with mandrel interior and polymer feet, suits for ND9 Vials
IPB250-P	29x5.7mm	250ul Plastic Insert with mandrel interior and polymer feet, suits for ND9 Vials
IA250	31x5mm	250ul Micro-Insert, clear glass, flat bottom, suits for ND8 Vials
IB300	31x6mm	300ul Micro-Insert, clear glass, flat bottom, suits for ND9 Vials
V18A	8.2x40mm	1ml shell vials, 1st hydrol glass, clear, 8mm PE-Plug, soft, without insertion