

Finetech®

NEW COMPACT SERIES

Oilless Vacuum Pump

VS300



Features

- Low vibration and low noise emission
- Oilless design, no pollution generated
- Compact and easy-to-carry design
- Rugged design and low maintenance
- Overload protection

VS400



Applications

- Laboratory instruments
- Medical instruments
- Use in pollution free environment
- Vacuum systems
- Vacuum filtration

Specifications

Model	Voltage	Open Flow (L/min)		Max. Vacuum (mmHg)	Power (W)	N. W. (kg)	HP
		50Hz	60Hz				
VS300	110v 60Hz or 220V 50Hz	19	24	660*	95	3.6	1/8
VS400	110v 60Hz or 220V 50Hz	32	40	660*	120	5.5	1/6

* Value \pm 6%

Dimension : 267x113x150 mm (VS300)

330x135x180 mm (VS400)

The Finetech logo consists of the word 'finetech' in a blue, lowercase, sans-serif font. To the right of the text is a stylized blue icon of a vacuum pump or a similar mechanical component.