

Qualitative And Quantitative Filter Paper

Centennial craft manufacturing

Made of 100% cotton pulp

By using only the finest raw materials in the production of its filter papers: cotton linters, high alpha cellulose pulp, borosilicate micro-glass and extremely pure micro-quartz fibers, we offer a wide range of filtration products to meet most every laboratory filtration requirement.



ITEMS	UNIT	Qualitative Filter Paper			Quantitative Filter Paper		
		Fast speed	Medium	Slow speed	Fast speed	Medium	Slow speed
		No.1	No.2	No.131	No.5A	No.5B	No.5C
Basis weight	GSM	80	80	80	80	80	80
Precipitate		Coarse and Gelatinous Precipitate	Crystalline	Fine Crystalline	Coarse and Gelatinous Precipitate	Crystalline	Fine Crystalline
Filter speed \leq	S	35	70	140	35	70	140
Ash Content \leq	%	0.15	0.15	0.15	0.01	0.01	0.01
Fe Content \leq	%	0.003	0.003	0.003	0.003	0.003	0.003
Water Extracts	PH	6.5-7.5	6.5-7.5	6.5-7.5	5-8	5-8	5-8
Pore Size	μm	20-25	15-20	10-15	20-25	15-20	10-15
Water Soluble Chloride \leq	%	0.02	0.02	0.02	0.02	0.02	0.02
Similar to Whatman brand		No.4	No.1	No.5	No.41	No.40	No.42 No.44
Size		60X60CM Diameter of Circular Filter Paper: 55.70.90.110.125.150.180MM We can supply other specifications according to the requirement of customers.			60X60CM Diameter of Circular Filter Paper: 55.70.90.110.125.150.180MM We can supply other specifications according to the requirement of customers.		



Finetech Research And Innovation Corp.

No.29, Anle St., Xiushui Township, Changhua County 504, Taiwan (R.O.C.)

TEL : 886-4-7687839 www.finetech-filter.com