


**Layer 1**  
Glass prefilter(GF)  
pore size from 1-10um



**Layer 2**  
Glass prefilter(GF)  
down to 1 um



**Layer 3**  
Membrane filter  
Nylon, PVDF, PTFE, etc.  
Filters down to 0.45um or 0.2um



# TriTech™

HPLC/Air/Vent Filters

## SYRINGE FILTERS

- HPLC Class
- Integrity Tested
- Easier to distinguish different membranes types(Metech)
- Multiple choices of Membranes
- Medical Grade Polypropylene Housing
- OEM Accepted

### Ordering Information

100pcs/pk

①	②	③	④	⑤
Membrane	Diameter	Sterile	Pore size	Membrane characteristics
PVDF PTFE Nylon MCE CA PES PP RC GF	2G          Metech 017=13mm 033=33mm	N(non-sterile) S(EO/Gamma)	022=0.22 um 045=0.45 um	I = hydrophilic O = hydrophobic

For example : PTFEG030N022O → 25mm, PTFE+GF, 0.22um

① ② ③ ④ ⑤