

10ml & 20ml Crimp-Top Headspace & Caps

- Applied to Shimadzu Varian, and other autosampler
- Caps are made of high quality aluminum



Specifications

Model No.	Size	Description
VH10A	22.5x46mm	10ml crimp headspace clear vial, flat bottom
VH10B	22.5x46mm	10ml crimp headspace amber vial, flat bottom
VH20A	22.5x75.5mm	20ml crimp headspace clear vial, flat bottom
VH20B	22.5x75.5mm	20ml crimp headspace amber vial, flat bottom
SC20B20A	Φ20mm	nature PTFE/nature silicone septa, 20mm crimp-top aluminum cap, 10mm centre hole
SBC20B20B	Φ20mm	nature PTFE/nature silicone speta, 20mm crimp-top blue magetic aluminum cap.10mm centre hole

100pcs/pk