

2ml Crimp-Top Vials & Caps

- Applied to Shimadzu Varian, and other autosampler
- Caps are made of high quality aluminum



Specifications

Model No.	Size	Description
V11A	11.6x32mm	2ml wide opening crimp-top vial,clear
V11B		2ml wide opening crimp-top vial with write-on spot,clear
V11C		2ml wide opening crimp-top vial,amber
V11D		2ml wide opening crimp-top vial with write-on spot,amber
SC11A11B	Φ12x32mm	white PTFE/red Silicone septa, 11mm crimp-top aluminum cap,5.5mm centre hole
SC11A11A	Φ11mm	white PTFE/red Silicone septa, 11mm crimp-top aluminum cap,5.5mm centre hole
SC11F11A-1		Non-slip Nature PTFE/White Silicone septa,11 mm crimp -top aluminum cap, 5.5 mm centre hole
SC11AA11A		pre-slit white PTFE/red Silicone septa, 11mm crimp-top aluminum cap,5.5mm centre hole

100pcs/pk